

In-gel Tryptic Digest for Protein ID by Mass Spectrometry

David Miyamoto 2/12/2002

This protocol is based on Shevchenko A, Wilm M, Vorm O, & Mann M. Mass spectrometric sequencing of proteins from silver-stained polyacrylamide gels. *Anal Chem* 1996, 68:850-8. I have used it with success on both Coomassie and silver-stained gel bands. The procedure includes reduction and acetamidation steps that may be skipped if desired. For heavily stained Coomassie bands, it is helpful to wash gel pieces for 1 hr in 100 mM NH_4HCO_3 prior to dehydrating with acetonitrile (step 2).

1. Excise band from Coomassie or silver stained gel. Cut gel band into 1 mm cubes using clean razor blade on a clean glass surface. Transfer to an Eppendorf tube.
2. Remove excess water with pipet. Add 25-35 μL acetonitrile to tube to cover gel pieces. Incubate 10 minutes at RT to dehydrate and shrink gel pieces.
3. Remove acetonitrile with pipet. Speed-vac to dryness for 10 minutes.
4. Swell gel particles in 150 μL 10 mM DTT in 100 mM NH_4HCO_3 . Incubate 1 hr at 56°C.
5. Cool to RT. Replace DTT solution with 150 μL 55 mM iodoacetamide in 100 mM NH_4HCO_3 . Incubate 45 minutes at RT in the dark with occasional vortexing.
6. Remove solution and wash gel pieces with 150 μL 100 mM NH_4HCO_3 . Incubate 10 minutes at RT.
7. Remove NH_4HCO_3 solution with pipet. Add 150 μL acetonitrile to dehydrate gel pieces. Incubate 10 minutes at RT.
8. Repeat wash steps 6 through 7. Remove acetonitrile and speed-vac to dryness for 10 minutes.
9. Place tubes in ice water bath and swell gel particles in 25-35 μL digestion buffer. Incubate 45 minutes in ice water bath. Digestion buffer consists of 12.5 ng/ μL trypsin (Promega sequence-grade

modified porcine trypsin, Cat. #V511A) in 50 mM NH_4HCO_3 . To make the digestion buffer, dissolve 20 μg Promega trypsin in 80 μL Promega trypsin buffer solution (50 mM acetic acid), and dilute with 50 mM NH_4HCO_3 to 12.5 ng/ μL .

10. Remove trypsin-containing buffer. Add 5-10 μL 50 mM NH_4HCO_3 without trypsin to keep pieces wet during cleavage. Incubate o/n at 37°C.

11. Spin 1' at 14,000 rpm to spin down gel pieces. Save supernatant in a separate PCR tube.

12. Add 20 μL 20 mM NH_4HCO_3 to cover gel pieces. Incubate 10 minutes at RT. Transfer supernatant to the PCR tube from step 11.

13. Add 25 μL 5% formic acid, 50% acetonitrile to the gel pieces. Incubate 20 minutes at RT.

14. Spin 1' at 14,000 rpm. Remove formic acid/acetonitrile solution and save in the same PCR tube from step 11.

15. Repeat formic acid extraction (steps 13 through 14) twice more.

16. Dry PCR tube in speed-vac to complete dryness. Store at -20°C until analysis.